

This sturdy, height adjustable back rest is designed to elevate and support the upper body at a range of different angles when seated in bed, allowing patients to comfortably read, eat meals or watch television.

Easily adjusts to five alternative heights, from a low lying angle to an upright sitting position.

Constructed from sturdy yet lightweight 25 mm aluminium frame with flexible mesh back rest for superior comfort,

airflow and support. The back rest is easy to clean and folds flat for convenient storage.

Can also be used with selected chairs for additional support or height adjustments.

CODE	BL9160
WIDTH	500 mm
HEIGHT	570 mm
DEPTH	390 mm
COLOUR	Blue
CAPACITY	110 kg
MATERIAL	Aluminium frame; breathable nylon mesh
WARRANTY	1year



