

# Etac Feed cutlery

Award-winning ergonomically designed aids for daily living



## Feed – Cutlery and Dose Mug

Feed feeding cutlery, child spoon and dose mug are excellent aids for both adults and children. The cutlery is designed to facilitate feeding, and the dose mug makes it easy to drink.

### Feed cutlery

This feeding cutlery is made of plastic with a flat and smooth shape that feels comfortable in the mouth. The shape makes it easy to feed a person that has difficulty opening their mouth.

The top of the handle is designed to provide a functional grip, and the length enables the person being fed to grip the cutlery.

Dimensions	Weight	Order Code	
220 mm	37 g	80403101	
200 mm	31 g	80403102	
220 mm	34 g	80403103	_
Care			
Dishwasher-saf	e		
	200 mm 220 mm	200 mm 31 g 220 mm 34 g	200 mm 31 g 80403102 220 mm 34 g 80403103 Care

#### Spoons, when children want to eat by themselves

When a child wants to feed themselves this spoon is indispensable.

The handle is slightly curved and oval with a support edge to prevent the hand from sliding down. The spoon is angle adjustable and comfortable in the mouth thanks to its flat shape. The strap is useful for children with poor gripping ability.

Description	Dimensions	Weight	Order Code
Feed Spoon, Adjustable Angled for Right Hand Angled for Left Hand	200 mm	40 g 40 g	80403301 80403302
Material	Care	40 g	00403302
White polyamide	Dishwasher-saf	e	

#### Dose mug

The dose mug has been designed to assist with drinking. The spout shape provides good contact with the lips and reduces the need to lean ones head back for secure pouring. This mug is transparent and comes with 50ml 1.7 oz and 100ml 3.4 oz markings.

Description	Height	Weight	Order Code	
Feed Dose Mug	70 mm	40 g	80404101	
Material	Care			
Polycarbonate	Dishwasher-	safe		







"Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances."

