



Etac Swift shower chair / stool

Designed to assist users and carers with every day routines.

 **etac**[®]
Creating Possibilities

Etac Swift shower stool/chair

Swift is a prize-winning series of shower stools/chairs designed with functionality and simplicity in mind. The basic shower stool can easily be transformed to a complete shower chair, depending on the user's needs.

Stable on any floor

Swift has a unique construction that adapts to uneven floors thanks to the flexible seat. The soft high quality ferrules adapt to the floor and ensures a good grip. All these features contribute to the users safety.

Changes with the circumstances

Swift is adaptable to different needs. It is possible to start using the shower stool, and add arm supports and back support as the users circumstances change. Swift is easy to adjust to the desired height with telescopic legs, no tools are needed. It can also be set one step lower at the front for a slightly forward tilted seat to aid standing up. Soft seats which are specially suited for showering are available should extra comfort be needed.



Easy to handle



Swift is easy to handle thanks to its light weight. Swift comes flatpacked and is easy to assemble without the need for any tools. It is also corrosion free.



Swift shower stool



Swift shower stool with arm supports.



Swift shower chair with arm supports



Easy to clean



Practical soap dish

With the soap dish your shampoo and soap is easy to reach. It also holds most shower nozzles.



Safe surface

The structure of the surface ensures that Swift doesn't become slippery when wet.



Swift is also available in a low version



Standard height

Description	Colour	Weight	Item no.	Dimensions
Swift Shower Stool	Green	3.1 kg 6.8 lbs	81701010 ✓	Dimensions Seat height: 42 – 57 cm 16½ – 22½" Seat height LOW: 34.5–42 cm 13½–16½" Total height: 80 – 95 cm 31½ – 37½" Total height LOW: 72.5–80 cm 28½–31½" Seat width: 54 21¼" Seat depth: 41cm 16¼" Total width: 54 cm 21¼" Total width, with arm supports: 54.5 cm 21" Total depth: 50 cm 19¾" Total depth LOW: 46.5 cm 18" Width between arm supports: 45 cm 17¾" Material: Polypropylene, Aluminium
Swift Shower Stool	Grey	3.1 kg 6.8 lbs	81701410 ✓	
Swift Shower Stool Low	Grey	2.9 kg 6.4 lbs	81701450 ✓	
Swift Shower Stool with Armrests	Green	3.9 kg 8.6 lbs	81701020 ✓	
Swift Shower Stool with Armrests	Grey	3.9 kg 8.6 lbs	81701420 ✓	
Swift Shower Chair with Armrests & Back	Green	4.7 kg 10.4 lbs	81701030 ✓	
Swift Shower Chair with Armrests & Back	Grey	4.7 kg 10.4 lbs	81701430 ✓	
Swift Shower Chair with Armrests & Back Low	Grey	4.5 kg 8.8 lbs	81701460 ✓	



Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.

Accessories



Back support	Item no.
Green	81701120
Grey	84005069



Back support pad	Item no.
Green	81701130
Grey	84005072

The back support pad is specially suited for showering since it doesn't become slippery when wet and gives an instantly warm feeling.
Material: Cellular polyethylene.



Soap dish	Item no.
Green	81701140
Grey	80209266



Seat pad	Item no.
Green	81701132 ✓
Grey	84005073 ✓

The seat pad is specially suited for showering since it doesn't become slippery when wet and gives an instantly warm feeling.
Material: Cellular polyethylene.




Arm supports pair	Item no.
Green	81701110
Grey	84005074



 **etac**®
Creating Possibilities

“Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances.”

 = Standard Stock Item