

# Etac Cloo toilet seat raiser

Designed to assist users and carers with every day routines.



# Etac Cloo toilet seat raiser

Cloo is a height adjustable toilet seat raiser which combines flexibility and simplicity. It has a discreet and soft design to blend in with your bathroom interior.

## Adjustable angle and height

Cloo can be adjusted in three different heights 6 cm, 10 cm or 14 cm. It can also be set with an anterior tilt which facilitates sitting down and standing up for users with reduced mobility in hip or knee joints.

### **Functional arm supports**

Cloo is available with or without detachable arm supports. The arm supports are foldable but can be fixed in a folded down position if required. Cloo can be used with only one arm support if space around the toilet is limited.

#### Easy to clean and mount

To aid cleaning Cloo has soft plain surfaces and the seat stays folded in the upright position. The seat can also be easily detached should this be preferred. Mounting is quick and tool-less.





With or without arm supports

Etac Cloo is available with or without arm supports.



#### Interior tilt

Etac Cloo can be set with an anterior tilt which facilitates sitting down and standing up for users with reduced mobility in hip or knee joints.



Easy to clean

The seat can be easily detached for easy cleaning.



Accessory



Description Item no. Arm supports, pair 80301222

**Material** Polypropylene, powder coated steel tubing



Description Item no. Heightening spacers, 80301221 14 cm 5½ Standard height adjustability: 6 and 10 cm 2¼" and 4" as standard. These are needed if 14 cm 5½" is required.

Material Polyamide

<b>Description</b> Cloo Toilet Seat Raiser	<b>Dimensions</b> Seat width 40 cm 15¾" Total width 60 cm 23½"	<b>Weight</b> 3.3 kg 7.28 lbs	Item no. 80301204
Cloo Toilet Seat Raiser with Armrests	Seat width 40 cm 15¾" Total width 60 cm 23½" Height of Armrests from seat 21 cm 8¼" Width between Armrests 49 cm 19¼"	<b>4.8 kg</b> 10.58 lbs	80301214

Material

Clean the product with a solvent-free cleaning agent (pH level 6–9), Seat: Polypropylene or with a 70% disinfectant solution. Can be decontaminated at max. Attachment: Aluminium 85°C / 185°F if so required.













"Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances."

