



Etac SupportBelt

Tools to help and support in mobilization and transfer.

 **etac**[®]

Creating Possibilities

Support belt with handles

Etac SupportBelt is used to offer support and assistance during mobilization and transfers. The diagonally placed handles allow the caregiver to maintain a more ergonomic grip.

SupportBelt is a useful tool when helping a person to stand up or sit down, to support walking or to assist during transfer between bed, wheelchair, chair, etc. SupportBelt has non-slip material on the inside and an adjustable belt buckle.



Etac SupportBelt has diagonally placed handles that allows the caregiver to maintain a more ergonomic grip.



Intended Use

- Support stand up / sit down
- Support walking training
- Support transfer between bed, wheelchair, chair, etc
- Transfer from wheelchair to toilet

Colour code - size labelling

	XS
	S
	M
	L
	XL


Description	Width	Length	Item no.	Specifications
SupportBelt XS, 4 handles, waist 65-105 cm	9 cm 3.5"	105 cm 41.3"	IM8010509	Specifications Outside: Nylon Inside: Nylon with PU-coating Padding: Polyetylen Handles: Polyester strap with polypropylen Belt Buckle: Acetal, nylon Washing Instructions Can be disinfected with LT-CIN isopropylalcohol 72 %). Do not use fabric softener or bleach.
SupportBelt S, 5 handles, waist 65-105 cm	13 cm 5.1"	105 cm 41.3"	IM8010513	
SupportBelt M, 6 handles, waist 90-130 cm	13 cm 5.1"	130 cm 51.2"	IM8013013	
SupportBelt L, 8 handles, waist 125-150 cm	13 cm 5.1"	150 cm 59.1"	IM8015013	
SupportBelt XL, 11 handles, waist 145-170 cm	15 cm 5.9"	170 cm 66.9"	IM8017015	





 **etac**[®]
Creating Possibilities

“Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances.”

 = Standard Stock Item