



Tasty – Glass, Beaker and Plate

Tasty is a range of lightweight products with a functional and discreet design to adapt to a standard laid table. They make eating easier and are all dishwasher-safe.



Tasty Glass

The glass has a thick neck and a stable foot to aid drinking.

The lightness of this glass, plus the ability of the bowl to rest against the top of the hand when lifting makes it easy to use.

Tasty Plate

The plate has a raised edge to facilitate picking up food with one hand only. The underside has a rubber strip attached so that the plate is positioned securely. The plate is also suitable in a feeding situation.

Tasty Beaker

Tasty beaker with its big handle and supporting edge permits various grips. The large base makes the beaker stable and the holder insulates against hot liquid when the beaker is lifted with a two hand grip. The supporting edge substitutes for an additional big handle, also saving space in the kitchen cabinet. A lid with a spout is included.

Description	Capacity	Weight	Material	Order Code	
Tasty Glass	3 dl	75 g	Polycarbonate	80404001	
Tasty Plate		160 g	White, heat-resistant melamine	80404005	
Tasty Beaker	3 dl	125 g	Polycarbonate & Polypropylene	80404042	

Care Dishwasher-safe. Not suitable for microwaves. Will retain its shiny finish if washed by hand





"Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances."

