



Etac PediTurn

Tools for seated transfers between bed, wheelchair, toilet, etc.

 **etac**[®]

Creating Possibilities

Etac PediTurn

Swivel cushions for the floor

Etac PediTurn is a swivel cushion that facilitates transfer between sitting positions, such as bed to wheelchair. It is used when the user cannot move their feet during the transfer.



Intended Use

- Support moving feet during seated transfer
- Transfer between bed and wheelchair
- Transfer between wheelchair, chair, shower chair, toilet, commode
- Transfer in and out of a car



Description	Diameter	Order Code
PediTurn	400 mm	IM99836GT



Specifications

Top: Toughtek
Inside: Nylon
Underside: Toughtek
Middle: Polyetylen

Washing Instructions


Can be wiped off with a damp cloth.





 **etac**[®]
Creating Possibilities

“Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances.”

 = Standard Stock Item