DJ**Med**[™] Finger Tip **Pulse Oximeter**CMS50D





- SpO2 value display
- Pulse rate value display, bar graph display
- Pulse waveform display
- Adjustable screen brightness
- Low-voltage indicator
- Auto power off function when not used
- Uses 2 AAA batteries
- Dimension: 57(L) × 31(W) × 32(H) mm
 Weight: About 50g (with the batteries)
- Display Mode: 0.96" Dual-color OLED
- Screen Resolution: 128*64
- SpO2 Range: 0%~100%, (1% resolution)
- Accuracy: 70%~100%: ±2%
- PR Measuring Range: 30bpm~250bpm
- Accuracy: ±2bpm or ±2%

SKU	Product	Price
1426	FINGER TIP PULSE OXIMETER CMS50D	\$114.95
8528	PULSE OXIMETER CASE	\$7.95







web djmed.com.au email info@djmed.com.au fax (02) 4744 2420 nsw (02) 8317 5097

qld (07) 3067 7073 vic (03) 8657 7678 sa (08) 7111 2981 wa (08) 6388 8026 To order call, email fax or shop online

DJMed™ Rechargeable Finger Tip **Pulse Oximeter**



- SpO2 value display
- Pulse rate value display, bar graph display
- Pulse waveform display
- Adjustable screen brightness
- Pulse rate sound indication
- Low-voltage indicator
- Auto power off function when not used
- Uses 3.7v rechargeable lithium battery
- Real-time data can be transmitted to computer via USB
- Dimension: 57(L) × 32(W) × 30 (H) mm Weight: About 50g (with lithium battery)
- Display Mode: 1.3" 65K Color OLED
- Screen Resolution: 128*96
- SpO2 Range: 0%~100%, (1% resolution)
- Accuracy: 70%~100%: ±2%
- PR Measuring Range: 30bpm~250bpm

Accuracy: ±2bpm or ±2%

SKU	Product	Price
1433	CMS50E RECHARGEABLE PULSE OXIMETER	\$249



ARTG 197923

Three easy ways to order

Shop online at DearJane.com.au			Fill out the order form and: Email: info@dearjane.com.au Or Fax: (02) 4744 2420			Call our wonderful, helpful staff (02) 8317 5097 Please have your credit card ready		
NAME:				NAME ON THE	CARD:			
ADDRESS:				CREDIT CARD	NUMBER:			
CITY:	STATE:	POSTCODE:		EXPIRATION:		CVV	:	(34 SIA BI) PE 88
PHONE: ()	EMAIL:			SIGNATURE:				
ITEM NUMBER	ITEM DESCRIPTION, SIZE OR COLOUR				QUANTITY	UNIT PRICE	TOTAL PRICE	
All prices are GST inclusive. Flat rate courier delivery to anywhere in Aus for orders under 15kg. Additional freight charge may apply for bulky go					· ·	PACKING & DELIV	ERY \$ 14.95	