



Etac Clean24 mobile shower commode

Designed to assist users and carers with every day routines.

 **etac**[®]

Creating Possibilities

Etac Clean 24"

Clean 24" is easy to use and propel.

Easy manoeuvring


The rear wheels on Clean 24" are positioned to maximize the manoeuvrability which makes the chair easy to propel without compromising with the security.

Easy transfers

Side transfers can be performed safely due to sufficient seat space in front of the wheels. Furthermore the brake handles are positioned just underneath the seat not to obstruct transfer but still being very easy to reach.

Improved access

The seat has a recess at the front and is fully open at the rear which enables unassisted intimate hygiene.



The brake handle is easy to reach without obstructing transfers.

Description

Clean 55 cm 21 3/4" with 24" rear wheels, green 80209216
Clean 55 cm 21 3/4" with 24" rear wheels, grey 80209276

Item no.



Materials

Chassis: Powder-painted steel tubing
Seat/back support: Polypropylene
Arm & Foot supports: Polypropylene
Castor brackets: Polyamide
Swivel castors, with ball bearings: polyamide/TPE
Rear wheels: Polyamide, mounted in ball bearings

Exterior dimensions:

Total width: 69 cm 27 1/2"
Width between arm supports: 43.5 cm 17 1/4"
Seat width: 48 cm 19"

Weight

20 kg 44.1 lbs



red dot



Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.

Accessories



Description
Clean Soft Seat Pad, green 80209236 ✓
Clean Soft Seat Pad, grey 80209260 ✓
Heat-reflecting soft seat with good anti-slip properties even when wet.
Recess width: 17 cm 6¾"
Material Cellular polyethylene



Description
Soap dish, Green 81701140
Soap dish, Grey 80209266
Attached to the side of the chair. Built-in holder for the shower handle.



Description
Cushion, Protect 80209245
Air-inflated ROHO cushion with excellent pressure distribution properties. The air pressure is adjustable.
Recess width: 21 cm 8¼"
Material: Neoprene



Description
Clean Cross Bar 80209246 ✓
Attached to the arm supports.
Material Aluminium, cellular polyethylene & polypropylene



Description
Clean Comfort pad 80209248 ✓
A comfortable, hygienic seat pad which makes the hole a little smaller.
Recess width: 18 cm 7"
Material Polyurethane



Description
Clean Safety belt 80209034 ✓
Attached to the rear of the back support.
Material Polyester



Description
Clean Comfort seat 80209227 ✓
Etac Clean's softest Comfort seat. The Comfort seat is contoured for added security and comfort. Intimate recesses at both front and back provide good accessibility.
The seat can easily be added and removed and can be cleaned at temperatures up to 85°C.
Material: polyurethane.



Description
Heightening spacers 80209243
Raise the seat by five centimetres. Do not fit Clean 24" rear wheels or Clean Height Adjustable.
Material: Aluminium



Description
Seat cover 80209247
When a complete seat is required.
Material: Polyurethane



Description
Rear wheel adapter 80209242
Only fits Clean with a seat height of 55 cm 21¾".
Material: Aluminium, polyamide



Description
Clean Soft Back Pad, green 80209238 ✓
Clean Soft Back Pad, grey 80209261 ✓
Heat-reflecting soft back pad.
Material Cellular polyethylene



Description
Elevating leg support, right 80209252
Elevating leg support, left 80209253
Attachments included. Only fits Clean height adjustable.
Material leg support: ABS plastic



Description
Clean Comfort, Comfort cover 80209228
A soft backrest with integrated armrest pads. The Comfort cover provides warm, soft support while the front of the armrest is left free to allow a secure grip when standing. The cover is cleaned using an alcohol wipe.
Material: polyethylene.



Description
Flip up foot supports, pair 80209250
Removable, adjustable height. Attachments included. Only fits Clean height adjustable.
Material foot support: Polypropylene



Description	Item no.
Clean Pan, green	80209230 ✓
Clean Pan, grey	80209254
Clean Lid <small>with handle</small> , green	80209232
Clean Lid <small>with handle</small> , grey	80209256
Clean Pan & Lid <small>with Handle</small> , green	80209231 ✓
Clean Pan & Lid <small>with Handle</small> , grey	80209255 ✓
Material Polyamide	



Description	Item no.
Attachment for low mounting of pan / bucket, green x 4	80209235 ✓
Attachment for low mounting of pan / bucket, grey x 4	80209259 ✓
Fitted to the pan/bucket holder.	
Material Polypropene	



Description	Item no.
Pan Holder, green	80209233 ✓
Pan Holder, grey	80209257 ✓
Fits both Clean pan and regular hospital pans with a thin rim.	



Description	Item no.
80mm non-lockable castors, green x 2	80209241
80mm non-lockable castors, grey x 2	80209264 ✓
Reduce seat height by 5 cm.	
The non-Lockable castors do not fit the Clean 24" rear wheels nor Clean Height Adjustable.	
Material polyamide/TPE	



Description	Item no.
Bucket Holder, Green	80209234
Bucket Holder, Grey	80209258
Fits bucket with a thick rim.	
Material Polypropene	



Description	Item no.
80mm lockable castors, green x 2	80209240 ✓
80mm lockable castors, grey x 2	80209263 ✓





 **etac**®
Creating Possibilities

"Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances."

 = Standard Stock Item