



Vac Maxi

Portable Suction Unit

The Vac series pump is APEX Medical's most recent generation of portable suction units and has been designed to combine function with aesthetic appeal and meet the needs of today's hygiene and infection control requirements. The Vac Maxi suction unit is a powerful, fast and effective aspirator which is ideal for use in high vacuum, high flow applications.

Vac Maxi

Portable Suction Unit

Features ▼

- Smooth hygienic finish, water resistant, easy to clean and maintain.
- 46 LPM powerful suction on demand.
- Optional accessories available

Single-patient use canister available in 800ml or 1800ml.

Filter and Tubing.



Integrated Canister Holder.



Portable design with ergonomic easy-to-carry handle.



Pressure regulator and easy-to-read analogue pressure gauge.



Detachable power cord allows you to quickly change to another cord to match local electrical outlet requirement.



Single-patient use canister available in 800ml or 1800ml.

Specification

Electrical Rating	AC 110 - 120V, 60Hz
	AC 220 - 240V, 50 / 60Hz
Power Consumption	100 W
Maximum Pressure	600 mmHg (adjustable)
Maximum Flow Rate	46 LPM
Net Weight	4.4 kg
Dimension	352 (L) x 206 (W) x 192 (H) mm
Sound Level	< 60 dB

DearJane Medical
orders@dearjane.com.au
DearJane.com.au
SKU: 4072

